

---

# BrainFlow Documentation

**Andrey Parfenov**

**Mar 30, 2020**



	<b>111</b>
.....	111
.....	111
.....	111
.....	111
.....	112
.....	112
.....	113
	<b>115</b>
.....	115
.....	115
.....	115
	<b>117</b>
	<b>119</b>
	<b>121</b>

BrainFlow is a library intended to obtain, parse and analyze EEG, EMG, ECG and other kinds of data from biosensors.

It provides a **uniform data acquisition API for all supported boards**, it means that you can switch boards without any changes in code and applications on top of BrainFlow are board agnostic. Also there is **powerful API to perform signal processing** which you can use even without BCI headset. Both of these two APIs are the same across bindings.



Table 1: Required inputs

Board	Board Id	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board	BrainFlow Board
Streaming Board	BoardIds.STREAMING_BOARD (-2)			multicast IP address	port	.	Board Id of master board	.
Synthetic Board	BoardIds.SYNTHETIC_BOARD (-1)			.	.	.	.	.
Cyton	BoardIds.CYTON_BOARD (0)	serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)	.	.	.	.	.	.
Ganglion	BoardIds.GANGLION_BOARD (1)	serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)	serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)	.	.	.	.	Timeout for device discovery(default 15sec)
Cyton Daisy	BoardIds.CYTON_DAISSY_BOARD (2)	serial port(COM3, /dev/ttyUSB0, /dev/cu.usbserial-xxxxxx...)	.	.	.	.	.	.
Ganglion Wifi	BoardIds.GANGLION_WIFI_BOARD (4)			Wifi Shield IP(default 192.168.4.1)	any local port which is free	.	.	.
Cyton Wifi	BoardIds.CYTON_WIFI_BOARD (5)			Wifi Shield IP(default 192.168.4.1)	any local port which is free	.	.	.
Cyton Daisy Wifi	BoardIds.CYTON_DAISSY_WIFI_BOARD (6)			Wifi Shield IP(default 192.168.4.1)	any local port which is free	.	.	.
BrainBit	BoardIds.BRAINBIT_BOARD (7)			.	.	.	Serial Number of BrainBit device	Timeout for device discovery(default 15sec)

## 1.1 Streaming Board

BrainFlow's boards can stream data to different destinations like file, socket and so on. This board acts like a consumer for data streamed from the main process.


**To use it in the first process you should call:**

```
# choose any valid multicast address(from "224.0.0.0" to "239.255.255.255") and port
start_stream (450000, 'streaming_board://225.1.1.1:6677')
```

- communication: None

## 1.3 OpenBCI

### 1.3.1 Cyton



[Cyton Getting Started Guide from OpenBCI.](#)

To choose this board in BoardShim constructor please specify:

- board\_id: 0
- serial\_port field of BrainFlowInputParams structure

Supported platforms:


- Windows >= 8.1
- Linux
- MacOS

**On MacOS there are two serial ports for each device: /dev/tty..... and /dev/cu..... You HAVE to specify /dev/cu.....**

Board Spec:

- num eeg(emg,...) channels: 8
- num acceleration channels: 3
- sampling rate: 250
- communication: serial port
- signal gain: 24

### 1.3.3 Cyton Daisy



CytonDaisy Getting Started Guide from OpenBCI.

To choose this board in BoardShim constructor please specify:

- board\_id: 2
- serial\_port field of BrainFlowInputParams structure

Supported platforms:

- Windows >= 8.1
- Linux
- MacOS


**On MacOS there are two serial ports for each device: /dev/tty..... and /dev/cu..... You HAVE to specify /dev/cu.....**

Board Spec:

- num eeg(emg,...) channels: 16
- num acceleration channels: 3
- sampling rate: 125
- communication: serial port
- signal gain: 24



### 1.3.4 Ganglion with WiFi Shield



[WiFi Shield Getting Started Guide from OpenBCI.](#)

[WiFi Shield Programming Guide from OpenBCI.](#)

To choose this board in BoardShim constructor please specify:

- board\_id: 4
- ip\_address field of BrainFlowInputParams structure should contain WiFi Shield Ip address(in direct mode its 192.168.4.1)
- ip\_port field of BrainFlowInputParams structure should be any local port which is free right now

Supported platforms:

- Windows >= 8.1
- Linux
- MacOS

Board Spec:

- num eeg(emg,...) channels: 4
- num acceleration channels: 3
- sampling rate: 1600
- communication: TCP socket to read data and HTTP to send commands

## 1.4 NeuroMD

### 1.4.1 BrainBit



BrainBit website.

To choose this board in BoardShim constructor please specify:

- board\_id: 7
- other\_info field of BrainFlowInputParams structure should contain Serial Number of BrainBit device
- optional: timeout field of BrainFlowInputParams structure, default is 15sec

Supported platforms:

- Windows  $\geq 10$
- MacOS

Board Spec:

- num eeg channels: 4
- num acceleration channels: None
- sampling rate: 250
- communication: Bluetooth Low Energy